

Resonance Ionization Mass Spectrometry

Reto Trappitsch




Brandeis
UNIVERSITY

February 28, 2022

RIMS — A Versatile Technique for Trace Element Analyses

- High sensitivity for small, atom-limited samples
- Minimal sample preparation
- Resonance ionization with tunable Ti:Sapphire lasers
- High spatial resolution
 - $\sim 1\ \mu\text{m}$ for laser desorption
 - $< 100\ \text{nm}$ for ion sputtering
- High useful yield
 - 38% for U analysis (Savina+ 2018)
 - $\sim 18\%$ for Ti analysis (Trappitsch+ 2018)
- Low backgrounds and high isobar suppression




A RIMS Table of Elements

A RIMS Periodic Table																	
H																	He
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	*	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	**															
			* La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
			** Ac	Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr


■ accessible by RIMS
■ published RIMS studies
■ published RIMS isotopic measurements

Savina and Trappitsch (2021)


An overview of Resonance Ionization Mass Spectrometry (RIMS)




An overview of Resonance Ionization Mass Spectrometry (RIMS)




An overview of Resonance Ionization Mass Spectrometry (RIMS)




An overview of Resonance Ionization Mass Spectrometry (RIMS)




An overview of Resonance Ionization Mass Spectrometry (RIMS)





An overview of Resonance Ionization Mass Spectrometry (RIMS)



An overview of Resonance Ionization Mass Spectrometry (RIMS)








Measurement Cycles repeat at 1 kHz


- 1 Desorption / Sputtering of sample
- 2 Ejection of secondary ions
- 3 Resonance ionization of photoions
- 4 Extraction
- 5 Mass / Charge separation and detection

Optional second ionization laser pulse allows for separation of isobars



Savina and Trappitsch (2021)

Sample Removal: Sputtering vs. Laser Desorption




- Sputtering with Ga ion beam
 - < 50 nm spatial resolution
 - Motionless blanking required
 - Trade off between high current or high spatial resolution
 - Duty cycle compared to SIMS: $\sim 10^{-4}$
- Desorption laser
 - Various wavelength possible to couple with different materials
 - Spot-size down to around $1\ \mu\text{m}$
 - Very low secondary ion backgrounds can be achieved

Sample Removal: Sputtering vs. Laser Desorption




EKSPLA 1064 nm
Desorption Laser



- Sputtering with Ga ion beam
 - < 50 nm spatial resolution
 - Motionless blanking required
 - Trade off between high current or high spatial resolution
 - Duty cycle compared to SIMS: $\sim 10^{-4}$
- Desorption laser
 - Various wavelength possible to couple with different materials
 - Spot-size down to around 1 µm
 - Very low secondary ion backgrounds can be achieved

Ionizing of Neutral Atoms: You only get One Chance!




Savina and Trappitsch (2021)

- Ionization laser beam size: ~ 1.5 mm diameter cylinder
- Laser intercepts cloud of neutrals above sample surface
- Neutrals that do not get ionized in first shot will be lost due to cloud expansion


Multi-step Laser Ionization of, e.g., Titanium (Trappitsch et. al, 2018)

- Resonance Ionization of Titanium requires three lasers
- Each ionization step is highly selective
- Ionization schemes need to be tested:
 - Spectroscopy of states above ionization potential
 - Saturation: Irradiance counts!
- Ti has low lying states
 - Understand population of these states
 - Scheme specific
 - Here: majority after sputtering in ground state




Multi-step Laser Ionization of, e.g., Titanium (Trappitsch et. al, 2018)

- Resonance Ionization of Titanium requires three lasers
- Each ionization step is highly selective
- Ionization schemes need to be tested:
 - Spectroscopy of states above ionization potential
 - Saturation: Irradiance counts!
- Ti has low lying states
 - Understand population of these states
 - Scheme specific
 - Here: majority after sputtering in ground state




Multi-step Laser Ionization of, e.g., Titanium (Trappitsch et. al, 2018)

- Resonance Ionization of Titanium requires three lasers
- Each ionization step is highly selective
- Ionization schemes need to be tested:
 - Spectroscopy of states above ionization potential
 - Saturation: Irradiance counts!
- Ti has low lying states
 - Understand population of these states
 - Scheme specific
 - Here: majority after sputtering in ground state





Multi-step Laser Ionization of, e.g., Titanium (Trappitsch et. al, 2018)

- Resonance Ionization of Titanium requires three lasers
- Each ionization step is highly selective
- Ionization schemes need to be tested:
 - Spectroscopy of states above ionization potential
 - Saturation: Irradiance counts!
- Ti has low lying states
 - Understand population of these states
 - Scheme specific
 - Here: majority after sputtering in ground state




Separation of m/q in Time-of-Flight (TOF) Mass Analyzer

- Time of Flight Mass Analyzer
 - ~ 3.5 m flight path
 - Grid-less reflectron to optimize transmission
 - Mass resolution $\frac{m}{\Delta m} > 1000$
- Difficulty: Map a photoion volume in time onto detector
- Lasers however take care of isobars




Ion Counting — Record every Arrival Time


- Ion counting detectors
 - Microchannel plate detectors (MCPs)
 - TOF Electron Multipliers
- Time-to-Digital Conversion: 80 ps time resolution
- Overall system dead-time: ~ 700 ns
- Reasonable count rates: $\sim 2,000$ cps




Simultaneous Measurements of Iron and Nickel



Simultaneous Measurements of Iron and Nickel




Multi-Element Analysis avoiding isobaric overlap



Full Disclosure — Limitations of RIMS

- Count rate limitations significantly limits the dynamic range
 - Narrowband lasers can be used in special cases to increase dynamic range
 - Example: $^{41}\text{Ca}/^{40}\text{Ca}$ analysis
- Ionization laser pulse width ~ 20 ns:
 \rightarrow Duty cycle $\sim 10^{-5}$
- Desorption laser coupling depends on material and wavelength
 - Choose the right wavelength and pulse width
- Sample material is removed as molecule
 - In-vacuo surface chemistry



Denk et al. (2006)

Full Disclosure — Limitations of RIMS

- Count rate limitations significantly limits the dynamic range
 - Narrowband lasers can be used in special cases to increase dynamic range
 - Example: $^{41}\text{Ca}/^{40}\text{Ca}$ analysis
- Ionization laser pulse width ~ 20 ns:
→ Duty cycle $\sim 10^{-5}$
- Desorption laser coupling depends on material and wavelength
 - Choose the right wavelength and pulse width
- Sample material is removed as molecule
 - In-vacuo surface chemistry



Full Disclosure — Limitations of RIMS

- Count rate limitations significantly limits the dynamic range
 - Narrowband lasers can be used in special cases to increase dynamic range
 - Example: $^{41}\text{Ca}/^{40}\text{Ca}$ analysis
- Ionization laser pulse width ~ 20 ns:
→ Duty cycle $\sim 10^{-5}$
- Desorption laser coupling depends on material and wavelength
 - Choose the right wavelength and pulse width
- Sample material is removed as molecule
 - In-vacuo surface chemistry



Hiden Analytical IG20

The Next Generation RIMS Instrument

- TOF Mass Spectrometer Optimization
 - Commercial TOF?
 - Optimized home-built TOF?
- Improved laser design and automation
- Ion Imaging
- Cryo-capability to handle biological samples
 - Trace isotopes in tissues, ...
 - Medical labeling with radioactive isotopes

**New capabilities, research areas, and
higher instrument up-time**




Kore Technology SurfaceSeer

The Next Generation RIMS Instrument

- TOF Mass Spectrometer Optimization
 - Commercial TOF?
 - Optimized home-built TOF?
- Improved laser design and automation
- Ion Imaging
- Cryo-capability to handle biological samples
 - Trace isotopes in tissues, ...
 - Medical labeling with radioactive isotopes

**New capabilities, research areas, and
higher instrument up-time**



Thank you!

Acknowledgement



THE UNIVERSITY OF
CHICAGO